



UNIVERSITY OF WISCONSIN-MADISON
UNIVERSIDAD NACIONAL AUTÓNOMA DE MÉXICO

AIUM 2018 ANNUAL CONVENTION
New York, NY, March 24-28, 2018

Quantitative Ultrasound Assessment of Neurotoxicity in the Young Rhesus Macaque Brain

Ivan M. Rosado-Mendez,^{1,2} Laura Castañeda,¹ Chrysanthy Ikonomidou,³
James A. Zagzebski,² and Timothy J. Hall²

¹*Instituto de Física, Universidad Nacional Autónoma de México, Mexico City, MEX*

²*Department of Medical Physics, University of Wisconsin, Madison, WI*

³*Department of Neurology, University of Wisconsin, Madison, WI*



Disclosures

The authors certify that they have NO financial interest in the subject matter or materials discussed in this presentation

Equipment loan and technical support from **SIEMENS Healthcare**

This work was supported by



GRANTS R01HD072077, R01HD083001-01A1
WNPRC PILOT AWARD, NIH P51OD011106



SEED FUNDING FROM INSTITUTO DE FISICA
UNAM PAPIIT IA104518 AND IN107916



Background

Every year, millions of children are exposed to anesthesia worldwide¹



Image from www.scientificamerican.com

FDA

WARNING

Prolonged and repeated use of general anesthetics in children younger than 3y/o and pregnant women can have neurotoxic effects²

Evidence based on *ex vivo* and *in vivo* studies in animal models³

Neurotoxic effects in humans have not been confirmed due to lack of methods to assess these effects non-invasively in the young brain

¹Andropolous and Greene, PMID: 28177852

²<https://www.fda.gov/downloads/Drugs/DrugSafety/UCM554644.pdf>


³McGowan and Davis, PMID: 18499584




Background

N-Methyl-D Aspartate (NMDA) antagonists and γ -Aminobutyric Acid A ($GABA_A$) agonists can induce apoptosis of neurons and oligodendrocytes

Apoptosis (programmed cell death) causes changes in tissue microstructure to which ultrasound is sensitive



Activated Caspase 3 (CASP3, an apoptotic enzyme) in neonatal Rhesus macaque (*Banbrick et al., PMID: 2222480*)



Oligodendrocytes stained for myelin basic protein after exposure of fetal Rhesus macaque to isoflurane (*Creeley et al., PMID: 24158051*)




Background

Ultrasound signals are sensitive to how tissue microstructure changes during apoptosis

This information can be obtained through the spectral (frequency-based) analysis of echo signals, i.e., Quantitative Ultrasound (QUS)⁴

QUS has been used *in vitro*, *ex vivo*, and *in vivo* to assess apoptosis and other types of cell death in various cancer types and tissues^{5,6}



Cell death evaluated through histology and QUS in breast tumor xenografts in mice (*Sadeghi et al., PMID: 23761215*)

⁴Czarnota, Kolios, et al. PMID: 10507779

⁵Vlad et al., PMID: 21057928

⁶Tadayyon et al. PMID:28710751



Objectives

LONG TERM GOAL

Investigate the use of Quantitative Ultrasound to detect the induction of apoptosis by anesthetics in the newborn brain

Current objectives:

1. Perform a preliminary assessment of the sensitivity of Quantitative Ultrasound (QUS) to apoptosis induced by sevoflurane (GABA_A agonist) in a non-human primate model of the newborn brain
2. Identify challenges inherent to the application of QUS techniques in the *in vivo* newborn brain




Methods

SIEMENS Acuson S3000

Phased array transducer @ 10MHz

- 14 2-7 day-old Rhesus macaques from the Wisconsin National Primate Research Center
- Five-hour long exposure to sevoflurane (GABA_A agonist)
- Ultrasound data acquisition at 0hrs and 6hrs after start of exposure
- Tissue harvested at 8hrs for activated caspase 3 immunohistochemistry




10-15 frames of “raw” radiofrequency (RF) signals from thalamus using Axius Direct User Research interfase⁷



QUANTITATIVE ULTRASOUND

DIFFUSE SCATTERING ANALYSIS


- Biomarkers for apoptosis -




1. Echo signal power spectrum

2. Attenuation compensation

3. Backscatter coefficient (BSC)



4. Assuming diffuse scattering,
ESS: BSC vs. frequency
AC: BSC magnitude



10-15 frames of "raw" radiofrequency (RF) signals from thalamus using Axis Direct User Research interface⁷



QUANTITATIVE ULTRASOUND

DIFFUSE SCATTERING ANALYSIS

- Biomarkers for apoptosis -

Effective Scatterer Size
(ESS)

Acoustic Concentration
(AC)

(Insana & Hall, PMID: 2299033)


REJECTION CRITERIA

Spectral signal/noise < 15 dB,
Disagreement between results from left and right thalamus

DATA ANALYSIS

6hrs. vs. 0hrs Parameter Change (P= ESS or AC)

$$\Delta P(\%) = \frac{P(6hrs) - P(0hrs)}{P(0hrs)} \times 100\%$$



10-15 frames of “raw” radiofrequency (RF) signals from thalamus using Axis Direct User Research interface⁷



QUANTITATIVE ULTRASOUND

DIFFUSE SCATTERING ANALYSIS

- Biomarkers for apoptosis -

Effective Scatterer Size
(ESS)

Acoustic Concentration
(AC)

(Insana & Hall, PMID: 2299033)

REJECTION CRITERIA

Spectral signal/noise < 15 dB,
Disagreement between results from left and right thalamus

DATA ANALYSIS

6hrs. vs. 0hrs Parameter Change (P= ESS or AC)

$$\Delta P(\%) = \frac{P(6hrs) - P(0hrs)}{P(0hrs)} \times 100\%$$

COHERENT SCATTERING ANALYSIS

- Validity of diffuse scattering assumption -


Detection of different scattering conditions
(Rosado-Mendez et al. PMID: 27046872)

Low
concentration

Diffuse
Scattering

Periodic
scatterers

Specular
reflectors






Results

Data rejection:

- Low Spectral signal/noise: 5 cases
- Disagreement left vs. right sides:
ESS: 3 cases, AC: 2 cases
- Remaining for analysis
ESS: 6 cases, AC: 7 cases

6hrs vs. 0hrs parameter change		
	ESS	AC
Mean	-5.4%	-2.6%
Range	-21.1% to 8.4%	-14.9% to 14.1%
Wilcoxon p	0.05	0.36
Significant reduction in:	5/6 subjects (83%)	4/7 subjects (57%)






Results

Data rejection:

- Low Spectral signal/noise: 5 cases
- Disagreement left vs. right sides:
ESS: 3 cases, AC: 2 cases
- Remaining for analysis
ESS: 6 cases, AC: 7 cases


6hrs vs. 0hrs parameter change		
	ESS	AC
Mean	-5.4%	-2.6%
Range	-21.1% to 8.4%	-14.9% to 14.1%
Wilcoxon p	0.05	0.36
Significant reduction in:	5/6 subjects (83%)	4/7 subjects (57%)






Coherent scattering leads to data paucity and increased variance^{8,9}

Total of measurements: 721



Average fraction of each subject's 10-15 frames that show a particular scattering condition



⁸Rubert et al., PMID 24726800

⁹Luchies et al. PMID: 22622974

*Agreement between QUS results from left and right thalamus



Conclusion

- The Effective Scatterer Size, a Quantitative Ultrasound parameter based on assuming diffuse scattering, shows promise as a biomarker for apoptosis induced by exposure to sevoflurane in the *in vivo* young Rhesus macaque brain
- The presence of coherent scattering conditions severely constrained our data. Current work is focusing on overcoming these limitations through data compounding and other strategies (Liu et al. PMID: 7521044, Herd et al. PMID: 22033132, Luchies et al. PMID: 22622974)



Acknowledgments

The authors thank:

Dr. Kevin Noguchi (Washington University at St. Louis)
for activated CASP3 analysis of harvested tissue

Kristin Crosno and staff at the Wisconsin National
Primate Research Center for coordinating experiments

Dr. Quinton Guerrero (University of Wisconsin)
for help in ultrasound data acquisition





Contact information

Ivan M. Rosado-Mendez, PhD
Research associate, Instituto de Física
Universidad Nacional Autónoma de México

irosado@fisica.unam.mx

